

## Push it



# Safety Syringe

All Sol-Care™ Safety Syringes are designed with an integrated safety feature which retracts the needle back into the barrel and allows you to break the plunger to prevent needle stick injury

**Sol-Care™ Safety Syringes with fixed needles** are ideal for insulin injections, TB testing and other immunizations

**Sol-Care™ Safety Syringe Allergy and Flu Trays** provide convenience for allergy and flu injections while reducing packaging waste

## Pull it



- *Permanently attached needle minimize medication waste. On average, 11 doses can be drawn from a 10 dose influenza vaccine multi-dose vial*

- *0.5ml, 1ml syringes facilitates small dose injections*

- *Potential reduction and prevention of disposal related injuries*

## Snap it





## SOL-CARE™ TB Safety Syringe with Fixed Needle



Ref	Size	Gauge	Length	Qty/Box	Qty/Case
100090IM	0.5mL	27	1/2"	100	1000
100018IM	1mL	25	5/8"	100	1000
100034IM	1mL	25	1"	100	1000
100070IM	1mL	26	3/8" IDB	100	1000
100019IM	1mL	27	1/2"	100	1000
100067IM	1mL	28	1/2"	100	1000
100082IM	1mL	30	1/2"	100	1000



## SOL-CARE™ Insulin Safety Syringe with Fixed Needle



Ref	Size	Gauge	Length	Qty/Box	Qty/Case
100002IM	0.5mL	29	1/2"	100	1000
100089IM	0.5mL	30	5/16"	100	1000
100071IM	1mL	28	1/2"	100	1000
100017IM	1mL	29	1/2"	100	1000
100081IM	1mL	30	1/2"	100	1000



## SOL-CARE™ Safety Trays



Ref	Size	Gauge	Length	Qty/Tray	Qty/Case
<b>Allergy Trays</b>					
100032IM	1mL	27	1/2"	25	1000
100069IM	1mL	26	3/8"	25	1000
<b>Flu Trays</b>					
100088IM	1mL	25	1"	25	1000