

Push it



Pull it



Snap it



Safety Syringe

All Sol-Care™ Safety Syringes are designed with an integrated safety feature which retracts the needle back into the barrel and allows you to break the plunger to prevent needle stick injury

Sol-Care™ Safety Syringes with fixed needles are ideal for insulin injections, TB testing and other immunizations

Sol-Care™ Safety Syringes with exchangeable needles are available in a variety of needle gauges and lengths

Sol-Care™ Safety Syringes without needles are available in four sizes

Sol-Care™ Safety Syringe Allergy and Flu Trays provide convenience for allergy and flu injections while reducing packaging waste

SOL-CARE™ TB Safety Syringe w/Fixed Needle

Ref	Size	Gauge	Length	Qty/Box	Qty/Case
100090IM	0.5mL	27	1/2"	100	1000
100018IM	1mL	25	5/8"	100	1000
100034IM	1mL	25	1"	100	1000
100070IM	1mL	26	3/8" IDB	100	1000
100019IM	1mL	27	1/2"	100	1000
100067IM	1mL	28	1/2"	100	1000
100082IM	1mL	30	1/2"	100	1000

