


 SOL-M™ Syringes & Needles

Sol-M syringes are available in a variety of sizes



Sol-M syringe and needle combinations deliver a comprehensive range of options

Sol-M hypodermic needles are available in a variety of gauges and lengths



Sol-M blunt fill needle and blunt fill needle with filter were developed specifically with medication preparation in mind


SOL-M[™] 1mL Syringe

Ref	Description	Size	Gauge	Length	Qty/Box	Qty/Case
180001PP	Syringe only LL tip	1mL	—	—	100	1000
180011	Syringe only slip tip	1mL	—	—	100	1000
181027	Tuberculin syringe	1mL	27	1/2"	100	800
181027S	Slip tip w/ exchangeable needle	1mL	27	1/2"	100	800
181025	Tuberculin syringe	1mL	25	5/8"	100	800
181025S	Slip tip w/ exchangeable needle	1mL	25	5/8"	100	800


SOL-M[™] 3mL Syringe

Ref	Description	Size	Gauge	Length	Qty/Box	Qty/Case
180003	Syringe only LL tip	3mL	—	—	100	800
1832558	Syringe w/ exchangeable needle	3mL	25	5/8"	100	800
1832510	Syringe w/ exchangeable needle	3mL	25	1"	100	800
1832310	Syringe w/ exchangeable needle	3mL	23	1"	100	800
1832210	Syringe w/ exchangeable needle	3mL	22	1"	100	800
1832215	Syringe w/ exchangeable needle	3mL	22	1 1/2"	100	800
1832110	Syringe w/ exchangeable needle	3mL	21	1"	100	800
1832115	Syringe w/ exchangeable needle	3mL	21	1 1/2"	100	800
1832010	Syringe w/ exchangeable needle	3mL	20	1"	100	800
1832015	Syringe w/ exchangeable needle	3mL	20	1 1/2"	100	800
1831815	Syringe w/ exchangeable needle	3mL	18	1 1/2"	100	800
1831815B	Syringe w/ blunt fill needle	3mL	18	1 1/2"	100	800


SOL-M[™] 5mL Syringe

Ref	Description	Size	Gauge	Length	Qty/Box	Qty/Case
180005	Syringe only LL tip	5mL	—	—	100	800
1852210	Syringe w/ exchangeable needle	5mL	22	1"	100	800
1852215	Syringe w/ exchangeable needle	5mL	22	1 1/2"	100	800
1852110	Syringe w/ exchangeable needle	5mL	21	1"	100	800
1852115	Syringe w/ exchangeable needle	5mL	21	1 1/2"	100	800
1852010	Syringe w/ exchangeable needle	5mL	20	1"	100	800
1852015	Syringe w/ exchangeable needle	5mL	20	1 1/2"	100	800
1851815B	Syringe with blunt fill needle	5mL	18	1 1/2"	100	800



SOL-M™ 10mL Syringe



Ref	Description	Size	Gauge	Length	Qty/Box	Qty/Case
180010	Syringe only LL tip	10mL	—	—	100	1200
1812210	Syringe w/ exchangeable needle	10mL	22	1"	100	1200
1812110	Syringe w/ exchangeable needle	10mL	21	1"	100	1200
1812115	Syringe w/ exchangeable needle	10mL	21	1 1/2"	100	1200
1812010	Syringe w/ exchangeable needle	10mL	20	1"	100	1200
1812015	Syringe w/ exchangeable needle	10mL	20	1 1/2"	100	1200
1811815B	Syringe w/ blunt fill needle	10mL	18	1 1/2"	100	1200



Ref	Description	Size	Gauge	Length	Qty/Box	Qty/Case
180020	Syringe only LL tip	20mL	—	—	50	600



SOL-M™ 30mL Syringe

Ref	Description	Size	Gauge	Length	Qty/Box	Qty/Case
180030	Syringe only LL tip	30mL	—	—	50	400



SOL-M™ 60mL Syringe

Ref	Description	Size	Gauge	Length	Qty/Box	Qty/Case
180060	Syringe only LL tip	60ml	—	—	30	240
180060CT	Syringe catheter tip	60ml	—	—	30	240

 **SOL-M[™] Blunt Fill Needle**



Ref	Size	Gauge	Length	Qty/Box	Qty/Case
110022	—	18	1 1/2"	100	1000

 **SOL-M[™] Blunt Fill Needle with 5 micron filter**



Ref	Size	Gauge	Length	Qty/Box	Qty/Case
110022F	—	18	1 1/2"	100	1000

 **SOL-M[™] Hypodermic Needle**



Ref	Size	Gauge	Length	Qty/Box	Qty/Case
113005	—	30	1/2"	100	1000
112905	—	29	1/2"	100	1000
112705	—	27	1/2"	100	1000
112506	—	25	5/8"	100	1000
112510	—	25	1"	100	1000
112310	—	23	1"	100	1000
112210	—	22	1"	100	1000
112215	—	22	1 1/2"	100	1000
112110	—	21	1"	100	1000
112115	—	21	1 1/2"	100	1000
112010	—	20	1"	100	1000
112015	—	20	1 1/2"	100	1000
111815	—	18	1 1/2"	100	1000