



## SOL-CARE™ Safety Allergy Trays



Ref	Size	Gauge	Length	Qty/Tray	Qty/Case
100032IM	1mL	27	1/2"	25	1000
100069IM	1mL	26	3/8"	25	1000



## SOL-CARE™ Safety Syringe with Fixed Needle




Ref	Size	Gauge	Length	Qty/Box	Qty/Case
100070IM	1mL	26	3/8" IDB	100	1000
100019IM	1mL	27	1/2"	100	1000
100067IM	1mL	28	1/2"	100	1000



## SOL-M™ 1 mL Allergy Tray

Ref	Description	Size	Gauge	Length	Qty/Tray	Qty/Case
181023T	Allergy mixing tray	1mL	23	1/2"	25	1000
181026IDBT	Allergy syringe tray	1mL	26	3/8" IDB	25	1000
181026T	Allergy syringe tray	1mL	26	1/2"	25	1000
181126IDBT	Allergy syringe tray	1mL	26	1/2" IDB	25	1000
181027T	Allergy syringe tray	1mL	27	1/2"	25	1000

Sol-Millennium offers a variety of allergy needle solutions. Choose between individual blister packs or trays; safety syringes or non-safety syringes; regular bevel needles or intradermal bevel needles.



- **Regular bevel needle** for intramuscular and subcutaneous injections. Sharp needle point allows for low penetration for, reducing discomfort to patient



- **Intradermal bevel needle (IDB)** designed for low angle insertion just below the epidermis. Typically used for skin testing