



SOL-M™ Luer Adapter



Discover our SOL-M™ Luer Adapter accessory

- An ideal accessory to use with a blood collection set that features a holder for evacuated tubes
- Manufactured with medical grade stainless steel, polypropylene, polyisoprene rubber
- Fitted with a rubber sleeve for ensuring that no blood flows out after successful vein penetration or between changing tubes - no blood spills!
- Perforated label not only serves to simplify identification, but also acts as a seal of product integrity
- Not made with natural rubber latex or DEHP
- 3 year lifespan

Reference number	Description	Box/Case (pcs)
LA21920	SOL-M™ Luer Adapter 20G	100/1000

For more information, please contact:

CORP-BC-112-0813

Sol-Millennium Medical Inc.
315 Shawnee North Drive, Suite 100,
Suwanee, GA 30024
T: 404-973-2200
e-mail: info@sol-m.com

sol-m.com

